

Determining the Absolute Configuration of the α -Carbon in JNJ-A (a GPR40 Superagonist) by NMR Spectroscopy

2019 Eastern Analytical Symposium & Exposition: NMR Spectroscopy Instrumentation and Application
11:10 AM, Nov. 19th, 2019


Ronghui Zhou, Ph.D.

rzhou11@its.jnj.com




Discovery Chemistry, Janssen R&D
Spring House, PA



Chirality




Handedness




Chiral Drug: the Thalidomide Tragedy in the 1960s



Thousands of children around the world were born with severe birth defects.



(*R*)-
thalidomide



(*S*)-
thalidomide

The *S* enantiomer causes birth defects, while the *R* enantiomer is effective against morning sickness.

Diabetes is a Global Health Crisis



DIABETES


DIABETES IS
ON THE RISE

A row of stylized human figures in blue. The figures vary in height and are arranged in a slightly curved line, representing a diverse group of people.

422 MILLION
adults have diabetes

3.7 MILLION
deaths due to diabetes
and high blood glucose


1.5 MILLION
deaths caused
by diabetes



THAT'S 1 PERSON IN 11




GPR40 Agonism is a Novel Mechanism for T2D




Defossa and Wagner (2014) *Bioorganic and Medicinal Chemistry Letters*, 24, 2291-3000

- A “nutrient sensor” GPCR on the pancreatic beta cell
- Highly expressed in the pancreas, intestine, also human brain
- GPR40 activation induces insulin secretion, in a glucose-dependent manner
- Expected to lower glucose in patients with T2DM with low risk of **hypoglycemia**
- GPR40 agonist has additional benefit of **weight loss**


GPR40 Agonists: Conformationally Constrained Fatty Acid Mimetics




TAK-875
Ph 3, Terminated
partial agonist




AM-837
Ph 1 Discontinued
partial agonist



AM-1638
full agonist




BMS-986118
full agonist





Merck, WO2016022448
full agonist

JNJ-A is a GPR40 Superagonist




EC₅₀ –Calcium Mobilization (nM)
(Mean ± SD)

	hGPR40 (n=)	rGPR40 (n=)
JNJ-A	0.25 ± 0.13 (6)	0.34 ± 0.18 (6)
JNJ-O	0.94 ± 0.37 (10)	4.1 ± 2.4 (10)



ACS Med. Chem. Lett. **2019**, 10, 1, 16-21

The Absolute Configuration of α -carbon in JNJ-A



How to determine the absolute configuration of α -carbon: **R** or **S**?

Absolute Configuration: Analytical Techniques



Single-crystal
X-ray Diffraction



Vibrational
Circular Dichroism




Chiral chromatography



NMR

Absolute Configuration: NMR Methods



the (*R*)/(*S*) notation
Cahn-Ingold-Prelog Rules

Chemical Shift

Coupling Constant

Nuclear Overhauser Effect

Computation


Machine Learning

Molecular Structure Feature

Single chiral center or remote chiral centers

Chiral Auxiliary

Chiral derivatizing agent: Mosher's acid
Chiral solvating/shift agent
Alignment medium



Molecular Structure Feature



Scalar Coupling Constant

Nuclear Overhauser Effect

Conformational Analysis

With neighboring predefined chiral center

Scalar Coupling Constant & the NOE Effect





Karplus Equation

$$^3J_{HH} = J_0 \cdot \cos^2 \Theta - K$$

$$J_0 = 14 \text{ (90-180°), } J_0 = 10 \text{ (0-90°), } K = 0$$

Bothner-By equation


$$^3J_{HH} = 7 - \cos \Theta + 5 \cdot \cos 2\Theta$$



Magritek Application Note 5: Glucose Anomers

The Absolute Configuration of the α -carbon in JNJ-A


Acyclic & highly flexible carbon-carbon single bond




What can we do to determine the absolute configuration of α -carbon: **R** or **S**?

Chiral Auxiliary: CDA, CSA, Alignment Medium
Conformational Analysis & Relative Stereochemistry

1D ^1H NMR Fully Assigned



H-F Coupling: Anti-conformer Predominates



Large coupling constant (~30 Hz)
between H₂₄ and F₃₀








Figure 3. Variation of H-F coupling with dihedral angle

Anal. Chem. 1965, 37, 3, 403-405




DFT Calculation: Anti-conformer Predominates




B3LYP-D3/6-31G** optimization of 20 conformers from OPLS3e conformational search
PCM water


MM/QM Calculation: Torsion/Energy Profiles




COOH




OPLS3e
GB/SA water




DFT
B3LYP-D3
PCM water
(rigid scan)




COO-







Why Anti-conformer: Steric Hindrance?




Relative size (A-value):



Steric Hindrance: Energy Barrier?








The Absolute Configuration of the α -carbon in JNJ-B




H instead of CH_3
How to determine the absolute
configuration of α -carbon: **R** or **S**?

Anti-conformer: beyond Steric Hindrance





DFT Calculation: Anti-conformer Predominates




B3LYP-D3/6-31G** optimization of 13 conformers from OPLS3e conformational search
PCM water

Anti-conformer: beyond Steric Hindrance




Further investigation:


1. What's the driving force for the anti-conformer in this case?
2. What's the role of carboxylic acid?




>>



>




The Nuclear Overhauser Effect: *R* or *S* at α -carbon





R-configuration




VS



S-configuration



R-configuration Confirmed by XRD



XRD Obtained on a Bruker D8 Platinum¹³⁵ CCD diffractometer equipped with Cu K_α radiation (λ = 1.5478) at UCSD Crystallography Lab.

Summary

1. The absolute configuration: NMR consistent with XRD
2. The steric hindrance around the α - & β -carbons locks the anti-conformation between the β -proton H₂₄ and fluorine
3. Key information: scalar coupling constant and NOE
4. MM/QM: energy profile and population distribution
5. NMR: convenient and reliable alternative to XRD

Acknowledgements

Computational Chemistry

Maxwell D. Cummings, Ph.D.
(mcummin1@its.jnj.com)

Analytical Chemistry

Gregory C. Leo, Ph.D.
A. Kerem Bingol, Ph.D.
John A. Masucci, Ph.D.
Gary W. Caldwell, Ph.D.
Zhikai Shi, Ph.D.

Medicinal Chemistry

Hui Huang, Ph.D.
(hhuang2@its.jnj.com)
Mark R. Player, Ph.D.
James C. Lanter, Ph.D.
Sanath K. Meegalla, Ph.D.
Michael P. Winters, Ph.D.

Discovery Chemistry, Discovery Sciences
Discovery, Product Development & Supply
Janssen Research & Development
Spring House, PA